



Automatic Autoclave TRA-602

www.safezex.com | info@safezex.com

Description

Automatic Autoclave TRA-602 is a microprocessor-controlled tabletop electronic device with fully automatic, pre-programmed sterilization cycles and an LCD screen that displays continuous cycle information. It has a 16 L capacity, a maximum temperature of 134 °C, and a pressure of 210 MPa. It has an electro-polished stainless steel chamber and door. A simple user interface is very easy to use.

Features:

- Tabletop autoclaves that are full auto and have pre-programmed sterilization cycles
- The chamber and entrance are made of stainless steel and have an electro-polish finish
- Control system with high precision for perfect sterilization results
- The built-in independent steam generator ensures that there is always steam accessible for faster cycles
- The user interface is very friendly and can operate very easily
- The LCD screen constantly provides data about the cycle
- The open-type water tank is very simple to fill and clean
- Microprocessor-controlled electronic autoclave that is fully automated
- It is a thermal autoclave with a non-vacuum working method
- Each operating status and program could be displayed as an icon on the LCD screen

Applications:

It is a thermal autoclave with a non-vacuum working method

Each operating status and program could be displayed as an icon on the LCD screen

Specifications :

Capacity	16 L
Dimension	445×400×690mm
Heating Power	2000 W
Maximum Pressure	210 MPa
Maximum Temperature	134 ?
Maximum Timer	0 to 55 min
Number Of Trays	3
Power Supply	220V 50Hz/110V60Hz
Sterilization Room Size	Ø230×360
Weight	45 kg



W 84th Terrace Lenexa KS 66214, USA
Email: info@safex.com | Website: www.safex.com