



## Automatic Sterilizer TRSA-603

[www.safezex.com](http://www.safezex.com) | [info@safezex.com](mailto:info@safezex.com)

## **Description**

Automatic Sterilizer TRSA-603 is a Computer Control Auto Recycle Sterilizer with a Stainless Steel Design. It has a capacity of 75 L and a maximum pressure range of 0.145 to 0.165Mpa. It delivers accurate results within a temperature range of 50°C to 126°C. Equipped with a quick-opening door framework with a hand wheel and a mechanical safety-interlocking device.

## **Features:**

It has a quick-open door structure with a hand wheel and a mechanical safety-interlocking device

A computer-controlled auto-recycled sterilization mechanism is built into the device

The LED screen displays the current state of operation

The temperature range is 50°C to 126°C

Having with the period of 0 to 9999min

Auto-discharging due to overpressure 0.145 to 0.165 mmHg

Sterilization is automatically terminated with a beep

Stainless steel is the only material used in its construction

## Specifications :

<b>Capacity</b>	75 L
<b>Dimension</b>	57×57×130 cm
<b>Heating Power</b>	4.5 KW
<b>Maximum Pressure</b>	0.145-0.165Mpa
<b>Maximum Temperature</b>	126 ?
<b>Maximum Timer</b>	9999min
<b>Power Supply</b>	380V 50Hz
<b>Sterilization Room Size</b>	?340×425mm
<b>Weight</b>	77 kg



W 84th Terrace Lenexa KS 66214, USA  
Email: [info@safzex.com](mailto:info@safzex.com) | Website: [www.safzex.com](http://www.safzex.com)