



Balance TRB-604

www.safezex.com | info@safezex.com

Description

Balance TRB-604 is a compact device with high stability; the sensor uses a unique locking device with ac/dc level, an LCD backlight display, and a weighing reaction speed range of 0 to 3000 g with precision of 10 mg. It is intended to be used in conjunction with a printer or computer.

Features:

It has ac/dc, level, LCD backlight display, fast weighing reaction time, and high stability

The sensor employs a one-of-a-kind locking mechanism

Counting function, percentage, automatic calibration, multiple calibration options, calibration weights

With unit conversion, stability, memory function, animal weighing, full range gross weight reduction, and so on

With an RS232 interface for connecting to a printer or computer

Specifications :

Accuracy	10 mg
Minimum Reading	0.01g
Pan Size	φ130 mm
Range	0 to 3000 g
Size	310×235×150 mm
Weight	1 Kg



W 84th Terrace Lenexa KS 66214, USA
Email: info@safezex.com | Website: www.safezex.com