



Centrifuge TRC-601

www.safezex.com | info@safezex.com

Description

Centrifuge TRC-601 is a small device that can accommodate a variety of sample tubes. It has a maximum speed of 4000 rpm and an RCF of 1790× g. It has a rotor capacity of 20 ml×6. It is very easy to use, comes with a switch to continue operating and set parameters, and is well suited for qualitative analysis.

Features:

It is having with a Scope of supply that includes the rotor

It can hold various sample tubes

It is having with a Small footprint

Very simple to operate

Having a Step speed adjustment

It Suits for qualitative analysis

Specifications :

Dimension	Ø250×220 mm
Maximum RCF	1790×g
Maximum RPM	4000 RPM, step speed adjustment
Rotor Capacity	20 ml×6
Timer	0 to 30min
Voltage	110V 60Hz /220V 50Hz
Weight	3.5 Kg



W 84th Terrace Lenexa KS 66214, USA
Email: info@safezex.com | Website: www.safezex.com