



Fume Hood TRFH-603

www.safezex.com | info@safezex.com

Description

Fume Hood TRFH-603 is a ventilation cabinet that combines water, electricity, gas. The window is equipped with an automatic balance type. It is having with a 1.8M fume hood and having with a ventilation capacity is 450 to 1600 m³/h. A stainless steel water tank with liquid antiseptic on the bottom. Internal air discharges at a variable speed of 0.3m to 0.6m/s by the ventilation cabinet. A low-speed exhaust fan is included.

Features:

Water, electricity, gas, and air conditioning are all in the cabinet

Inside, there is a multifunctional socket

A stainless steel water tank with liquid antiseptic is installed on the bottom

It is secure and dependable

The ventilation capacity is 450 to 1600 m³/h, and the air circulation effectiveness is high

It has a sliding glass door that can be opened and closed

The flow plate makes bad gases at different heights discharged from different parts

Using a water channel to solve a problem

Comes with a 1.8 M fume cupboard

It has a cabinet for chemical fumes

It is an all-stainless steel fume cupboard

Specifications :

Dimension Of Working Area	1630×700×600 mm
Discharging Speed	within 0.3 m to 0.6 m/s adjustable
Fume Hood	1.8M
Inlet Pressure	>0.5 Pa
Overall Dimension W×D×H	1800×760×1950 mm
Power	800 W
Power Supply	AC, 220 V/50 Hz
Weight	200 Kg



W 84th Terrace Lenexa KS 66214, USA
Email: info@safezex.com | Website: www.safezex.com