



Digital Heating Water Bath TRHB-602

www.safezex.com | info@safezex.com

Description

Digital Heating Water Bath TRHB-602 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 6 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

It is equipped with a microprocessor control

mechanism with a stainless steel interior

Designed with a digital display inner temperature

It is covered with an aluminum ring

It is a thermostatic device

Specifications :

Capacity	6 L
Dimension	37x23x24 cm
Power	600 W
Temperature Fluctuation	±0.5 °
Temperature Homogeneity	±1 °
Temperature Range	(RT+5) -300 °
Voltage	110V 60Hz /220V 50 Hz
Weight	4.5 kg
Working Size	30x16x12 cm



W 84th Terrace Lenexa KS 66214, USA
Email: info@safex.com | Website: www.safex.com