



## Digital Heating Water Bath TRHB-604

[www.safezex.com](http://www.safezex.com) | [info@safezex.com](mailto:info@safezex.com)

## **Description**

Digital Heating Water Bath TRHB-604 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 17 L. It is designed with a stainless steel interior with an aluminum ring cover.

## **Features:**

It is equipped with a microprocessor control

mechanism with a stainless steel interior

Designed with a digital display inner temperature

It is covered with an aluminum ring

It is a thermostatic device

## Specifications :

<b>Capacity</b>	17 L
<b>Dimension</b>	53x37x24 cm
<b>Power</b>	1500 W
<b>Temperature Fluctuation</b>	$\pm 1.5$ ?
<b>Temperature Homogeneity</b>	$\pm 1$ ?
<b>Temperature Range</b>	(RT+5) -300 ?
<b>Voltage</b>	110V 60Hz /220V 50 Hz
<b>Weight</b>	7.5 kg
<b>Working Size</b>	46x30x12 cm



W 84th Terrace Lenexa KS 66214, USA  
Email: [info@safzex.com](mailto:info@safzex.com) | Website: [www.safzex.com](http://www.safzex.com)