



Digital Heating Water Bath TRHB-605

www.safezex.com | info@safezex.com

Description

Digital Heating Water Bath TRHB-605 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 22 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

It is equipped with a microprocessor control

mechanism with a stainless steel interior

Designed with a digital display inner temperature

It is covered with an aluminum ring

It is a thermostatic device

Specifications :

Capacity	22 L
Dimension	67x37x24 cm
Power	2000 W
Temperature Fluctuation	±1.5 °
Temperature Homogeneity	±1 °
Temperature Range	(RT+5) -300 °
Voltage	110V 60Hz /220V 50 Hz
Weight	9 Kg
Working Size	60x30x12 cm



W 84th Terrace Lenexa KS 66214, USA
Email: info@safex.com | Website: www.safex.com