



Mini Vertical Electrophoresis cell TRMV-601

www.safezex.com | info@safezex.com

Description

Mini Vertical Electrophoresis cell TRMV-601 has a capacity of 1-2 pieces and it is used for pre-casting and hand-casting gel. There are no more than four pieces of gel available at any given time. It has a multifunctional gel scraper built-in. The superior comb does not slow gel polymerization. The gel is separated from the air by a built-in kit. Ensures a high level of leak resistance. Supplied with an 8373mm gel size.

Features:

The filler strip is stable to the glass to avoid leakage of gel

A built-in kit separates the gel from the air

The gel-making kit safeguards accurate aligning

The distinctive comb prepares not to stop gel polymerization

Excellent leak proof ness

The thickness and well-marked comb are understandable

Equipped with a multi-functional gel scraper

Specifications :

Comb Specification	1.0mm: 10 and 15 wells
	0.75mm: 10 and 15wells
	1.5mm: 10 and 15 wells
Dimension (L ×W×H)	180×130×160mm
Gel Dimension	83×73mm
Quantity	1 to 2 piece(s)
Short glass Plate	100×73mm
Spacer glass plate	100×83mm
Thickness	1.0mm: 10, 15 wells (standard), 0.75mm: 10 & 15 wells, 1.5mm: 10 & 15 wells (optional)



W 84th Terrace Lenexa KS 66214, USA
Email: info@safex.com | Website: www.safex.com