



## **N Vertical Autoclave TRVA-601**

[www.safezex.com](http://www.safezex.com) | [info@safezex.com](mailto:info@safezex.com)

## Description

N Vertical Autoclave TRVA-601 is made entirely of stainless steel, is simple to use, safe, and dependable, and comes with two stainless steel sterilizing baskets. Presented with a 35 L chamber volume, the working temperature was 134 °C, and the working pressure was 0.22 MPa. If no beeping reminder, the machine will shut down automatically. A digital display of the working status, and a touch key panel.

## Features:

Design made completely of impeccable high-combination steel

Computerized show of working status, a dash of key

After cleansing, the cool air and steam are naturally released

After cleansing, the machine will naturally close down with a blare to remind you

Outfitted with two treated steel disinfecting containers

It is not tough to work, protected and reliable

## Specifications :

<b>Adjustment Of Temperature</b>	105 to134 °C
<b>Chamber Volume</b>	35 L
<b>G.W</b>	56 Kg
<b>Heat Average</b>	?±1 °C
<b>Maximum Working Pressure</b>	0.23 MPa
<b>N.W</b>	42 Kg
<b>Overall Dimension</b>	480×460×850 mm
<b>Power</b>	AC220V.50Hz
<b>Power Consumption</b>	2.5Kw
<b>Timer</b>	0 to 99 min or 0 to 99 hour 59 min
<b>Transport Dimension</b>	570×550×970 mm
<b>Working Pressure</b>	0.22 MPa
<b>Working Temperature</b>	134 °C



W 84th Terrace Lenexa KS 66214, USA  
Email: [info@safezex.com](mailto:info@safezex.com) | Website: [www.safezex.com](http://www.safezex.com)