



NC Linear Shaker TRLS-602

www.safezex.com | info@safezex.com

Description

NC Linear Shaker TRLS-602 has dual LCD screens for parameter setting and viewing. It can withstand a maximum load weight and operates in a timed and constantly functioning mode. The dual LCD screens separately display time and speed, with a speed range of 100 to 350 rpm and a time range of 1 minute to 19 hours 59 minutes. Present with the brushless DC motor is maintenance-free.

Features:

There is a maximum load-bearing weight of 7.5kg allowed

The period is from 1min 19 hrs. 59 min

The time and speed are displayed separately on dual LCD screens

There are numerous inserts and fittings available

Designed with a brushless DC motor requires no maintenance

Specifications :

Ambient Temperature	5 to 40°C
Amplitude	10mm
Humidity	0.8
Maximum Load	7.5kg
Motor Type	DC motor
Speed Display	LCD
Speed Range	100 to 350rpm
Time setting range	0 to 19h 59min
Voltage	100 to 240V, 50/60Hz



W 84th Terrace Lenexa KS 66214, USA
Email: info@safexex.com | Website: www.safexex.com