



## Portable Sterilizer TRPS-601

[www.safezex.com](http://www.safezex.com) | [info@safezex.com](mailto:info@safezex.com)

## Description

Portable Sterilizer TRPS-601 is a stainless steel device with an electric heating mode and a Dual scale numerical indication pressure gauge. It is having with a movable interior that is convenient for operation and with a capacity of 18 L and disinfecting Pressure of 0.14-0.165Mpa. It is having with the self-inflating type seal.

## Features:

It is made of stainless steel

It has an overpressure auto-discharging system. 0.145 to 0.165MPa

The maximum working temperature is 126 °C

Numeric indication with a dual magnitude Seal that self-inflates

There are two modes: electric heating and electric/coal heating

The movable interior makes operation easier

## Specifications :

<b>Capacity</b>	18 L
<b>Dimension</b>	39×39×54 cm
<b>Disinfecting Pressure</b>	0.14-0.165 MPa
<b>Heating Mode</b>	Electric-heating
<b>Maximum Working Temperature</b>	126 ?
<b>Power</b>	2 KW
<b>Voltage</b>	110 V 60 Hz /220 V 50 Hz
<b>Weight</b>	14 kg



W 84th Terrace Lenexa KS 66214, USA  
Email: [info@safezex.com](mailto:info@safezex.com) | Website: [www.safezex.com](http://www.safezex.com)