



## Ultrasonic Cleaner TRUC-603

[www.safezex.com](http://www.safezex.com) | [info@safezex.com](mailto:info@safezex.com)

## **Description**

Ultrasonic Cleaner TRUC-603 has an ultrasonic vibration frequency of interaction composed of high-quality SUS304 stainless steel stampings. It has a 6 L capacity, a heating temperature range of 20 to 80°C, and a digital timer with a setting range of 0 to 99 minutes. It does not have a drainage ball valve.

## **Features:**

The intensity of ultrasonic vibration, thoroughly solve the ultrasonic cleaning

Its work-performed tank is made of high-quality SUS304 stainless steel molding

The bladder is made of high-quality steel material, and the cladding is made of aluminum

Electronic timer with a period range of 0 to 99 minutes

The drainage ball valves are installed in drainage systems 6L and higher models

## Specifications :

|                             |             |
|-----------------------------|-------------|
| <b>1 Unit Packing</b>       | 43×28×42 cm |
| <b>1 Unit Weight</b>        | 6 Kg        |
| <b>Dimensions</b>           | 33×18×31 Cm |
| <b>Drainage</b>             | With        |
| <b>Heating Power</b>        | 300W        |
| <b>Heating Temperature</b>  | 20 To 80?   |
| <b>Overall Packing</b>      | 86×46×82 Cm |
| <b>Rinse Tank Size</b>      | 30×15×15 cm |
| <b>Timer</b>                | 0 To 99min  |
| <b>Ultrasonic Frequency</b> | 28/40khz    |
| <b>Ultrasonic Power</b>     | 180W        |
| <b>Volume</b>               | 6 L         |
| <b>Weight</b>               | 38 Kg       |



W 84th Terrace Lenexa KS 66214, USA  
Email: [info@safex.com](mailto:info@safex.com) | Website: [www.safex.com](http://www.safex.com)