



Water Distiller TRWD-601

www.safezex.com | info@safezex.com

Description

Water Distiller TRWD-601 is corrosion and aging resistant and has an intelligent automatic control. It is equipped with a stainless steel bowl with an electric heater. It condenses during the steaming process for the distilled water by electric heating. It has an output water level of 5 liters per hour.

Features:

It is corrosion-resistant and age resistant

Having a smart automatic control

Electric heating causes steam condensation into distilled water

It is outfitted with distillation generating equipment, which includes a condenser and other components

It consists of a stain-free steel electric double distiller

Specifications :

Output	5 L/ hour
Package	300×300×630 mm
Power	AC220 V/ 4.5 KW
Weight	7 kg



W 84th Terrace Lenexa KS 66214, USA
Email: info@safezex.com | Website: www.safezex.com