



Water Distiller TRWD-602

www.safezex.com | info@safezex.com

Description

Water Distiller TRWD-602 is corrosion and aging resistant and has an intelligent automatic control. It is equipped with a stainless steel bowl with an electric heater. It condenses during the steaming process for the distilled water by electric heating. It has an output water level of 10 liters per hour.

Features:

It is corrosion-resistant and age resistant

Having a smart automatic control

Electric heating causes steam condensation into distilled water

It is outfitted with distillation generating equipment, which includes a condenser and other components

It consists of a stain-free steel electric double distiller

Specifications :

Output	10 L/hour
Package	340×340×635 mm
Power	AC380V/ 7.5 KW
Weight	8 kg



W 84th Terrace Lenexa KS 66214, USA
Email: info@safezex.com | Website: www.safezex.com