



Water Steam Distiller TRWD-606

www.safezex.com | info@safezex.com

Description

Water Steam Distiller TRWD-606 is a stainless steel water tower system that is resistant to corrosion and aging. Distillation of power plants, condensers, evaporators, electric heating pipes, and controls. It has an automatic and intelligent control system. It is equipped with a control panel and has a button mode. Equipped with a water level window that can display the remaining water level. The heating component is composed of an immersion heating tube and is durable.

Features:

It is a water tower mechanism made from steel

It is corrosion-resistant and age resistant

It has an automatic and intelligent control system

designed with a water level window that displays the residual water level

The heating component is made up of immersion heating tubes that are long-lasting

Specifications :

Output	20L/hour
Package	740×400×680mm
Power	AC380V/13.5 KW
Weight	17 Kg



W 84th Terrace Lenexa KS 66214, USA
Email: info@safezex.com | Website: www.safezex.com